

Radiation Dosimetry and Organ Doses from Imaging

May 2011

Choonsik Lee
Radiation Epidemiology Branch
DCEG/NCI/NIH

Contents

- Introduction
 - Radiation imaging
 - Radiation dosimetry in imaging
 - Organ dose estimation
- Dosimetry in major imaging modalities
 - Radiography
 - Mammography
 - Fluoroscopy
 - Computed tomography


Radiation Imaging

X-ray

- Discovery
 - Discovered by Wilhelm Rontgen (1895)
("X" is indicating "unknown")
- Findings
 - Travels in straight lines
 - Make shadows of absorbing material on photosensitive paper!




X-ray generation




Different imaging modalities



Radiography



Mammography



**Nuclear
Medicine
Imaging***



Computed Tomography



Interventional fluoroscopy



*Use different mechanism from other imaging modalities


Current status of procedures (US 2006)*

Modality	Number of procedures	%
Radiography	293 million	74
CT	67 million	17
Nuclear Medicine	18 million	5
Interventional Fluoroscopy	17 million	4
Radiotherapy	1 million patients	NA

* Mettler et al. Radiology (2009)


Changes in U.S. medical radiation exposure*

1980



Total 3.6 mSv (effective dose) per capita

2006



Total 6.2 mSv (effective dose) per capita

Typical organ doses for patient*

Study type	Relevant organ	Organ dose (mGy)
Dental radiography	Brain	0.005
PA chest radiography	Lung	0.01
Lateral chest radiography	Lung	0.15
Screening mammography	Breast	3
Adult abdominal CT	Stomach	10
Barium enema	Colon	15
Neonatal abdominal CT	Stomach	20

* Brenner et al. NEJM (2007)

Typical organ doses for patient*

Type of examination	Effective dose (mSv) ^f
Radiography (single radiograph) ^a	
Skull AP or PA	0.015 (1)
Chest PA	0.013 (1)
T-spine AP	0.27 (20)
L-spine AP	0.44 (30)
Abdomen AP	0.46 (35)
Pelvis AP	0.48 (35)
Mammography (4 views) ^b	
Screening	0.2 (15)

* Linet et al. PR (2009)

Typical organ doses for patient*

Type of examination	Effective dose (mSv) ^f
Dental radiography ^c	
Intra oral	0.013 (1)
Panoramic	0.012 (1)
Diagnostic fluoroscopy procedures	
Barium swallow ^a	1 (70)
Barium meal ^a	2 (150)
Barium enema ^a	5 (350)
Angiography—cerebral ^c	2 (150)
Angiography—cardiac ^c	7 (500)

* Linet et al. PR (2009)

Typical organ doses for patient*

Type of examination	Effective dose (mSv) ^f
Computed tomography ^d	
Head	2 (150)
Chest	10 (750)
Abdomen	10 (750)
Pelvis	7 (500)
Abdomen/pelvis	15 (1,100)
C-spine	5 (400)
T-spine	8 (550)
L-spine	7 (500)

* Linet et al. PR (2009)

Typical organ doses for patient*


Type of examination	Effective dose (mSv) ^f
Diagnostic nuclear medicine ^e	
Bone (^{99m} Tc-phosphate)	3 (200)
Heart (²⁰¹ Tl thallous chloride)	13 (950)
Lung (^{99m} Tc-MAA)	0.9 (70)
Tumor-PET(¹⁸ F-FDG)	7 (500)
Kidney (^{99m} Tc-MAG3)	0.6 (40)
Thyroid (^{99m} Tc-Pertechnetate)	0.9 (70)

* Linet et al. PR (2009)

Radiation dosimetry in imaging

What is Dosimetry?

- Definition: determination of radiation dose resulting from the interaction of ionizing radiation with matter




The diagram illustrates radiation exposure pathways for three scenarios:

- Nuclear medicine:** A yellow silhouette of a patient is shown with a radiation symbol on the chest. Red arrows indicate exposure to the head, neck, and torso. A red arrow points from the patient to a gray silhouette of a radiologist, indicating exposure to the radiologist.
- Radiography, Fluoroscopy, Mammography, CT:** A radiation symbol is shown on the left. A red arrow points from it to a gray silhouette of a patient, indicating exposure to the patient. A red arrow points from the patient to a gray silhouette of a radiologist, indicating exposure to the radiologist.
- Family Members (Nuclear medicine):** A yellow silhouette of a patient is shown with a radiation symbol on the chest. Red arrows indicate exposure to the head, neck, and torso. A red arrow points from the patient to a gray silhouette of a family member, indicating exposure to the family member.

Dosimetric quantities*

Kinetic energy deposited in matter



Kinetic energy released in matter

Radiation weighting factor*

Table 2. Recommended radiation weighting factors.

Radiation type	Radiation weighting factor, w_R
Photons	1
Electrons ^a and muons	1
Protons and charged pions	2
Alpha particles, fission fragments, heavy ions	20
Neutrons	A continuous function of neutron energy (see Fig. 1 and Eq. 4.3)

$$w_R = \begin{cases} 2.5 + 18.2 e^{-[\ln(E_n)]^2/6}, & E_n < 1 \text{ MeV} \\ 5.0 + 17.0 e^{-[\ln(2E_n)]^2/6}, & 1 \text{ MeV} \leq E_n \leq 50 \text{ MeV} \\ 2.5 + 3.25 e^{-[\ln(0.04E_n)]^2/6}, & E_n > 50 \text{ MeV} \end{cases}$$

Tissue weighting factor*


Table 3. Recommended tissue weighting factors.

Tissue	w_T	$\sum w_T$
Bone-marrow (red), Colon, Lung, Stomach, Breast, Remainder tissues*	0.12	0.72
Gonads	0.08	0.08
Bladder, Oesophagus, Liver, Thyroid	0.04	0.16
Bone surface, Brain, Salivary glands, Skin	0.01	0.04
	Total	1.00

* Remainder tissues: Adrenals, Extrathoracic (ET) region, Gall bladder, Heart, Kidneys, Lymphatic nodes, Muscle, Oral mucosa, Pancreas, Prostate (♂), Small intestine, Spleen, Thymus, Uterus/cervix (♀).


Dosimetric quantities*

Kinetic energy deposited in matter



Kinetic energy released in matter





Organ dose estimation for medically-exposed patients



- Controlled
- Relatively well documented


How to obtain the individualized organ dose?

(1) Measurement



- Expensive
- Substantial man-hour
- Not individualized

(2) Calculation



- 30+ organ doses
- Bone marrow dose
- Highly individualized

- Cost-effective
- Fewer man-hour
- More flexible

(3) Conversion Factor



Measurement




Conversion
Factor

Organ doses




Derived from computer simulations

Dosimetry in major imaging modalities



Radiography



Factors affecting dose in radiography

- Beam energy
 - Primarily depends on the tube potential (kVp) and filtration
 - Higher energy beam is more penetrating to reach image receptor
 - Lower tube current or shorter imaging time
 - Reduce the dose to the patient
- Filtrations
 - Total filtration = Inherent filtration + Added filtration
 - Remove low-energy x-ray which can be absorbed by the patient
- Collimation
 - Limit the exposed area in the patient
 - Reduce the scattered radiation and increase image contrast

Factors affecting dose in radiography

- Grids
 - Reduce the scattered radiation contribution to improve image contrast
 - Also absorb a portion of non-scattered radiation
 - Cause increase current and time giving more doses to the patient
- Patient size
 - Need more radiation to get an acceptable image for thicker patient
 - Technique charts displaying suggested technique factors for different exams and patient thicknesses will be helpful

Organ dose estimation: Conversion factors

- “Handbook of selected tissue doses for projections common in diagnostic radiology” (Rosenstein, FDA89-8031, 1988)
 - Developed from adult male and female computational phantoms coupled with Monte Carlo transport technique
 - Provide organ doses per unit exposure (measurable) for comprehensive technique factors

Organ dose estimation: Conversion factors

TABLE 24. PA CHEST - SID: 72" (183 cm); FIELD SIZE at FILM: 14" X 17" (35.6 cm X 43.2 cm)

MALE HVL (mm Al) →	TISSUE DOSES (mrad) and CANCER DETRIMENT INDEX for 1 R EXPOSURE at SKIN ENTRANCE (FREE-IN-AIR) ^{a, b}											
	1.0	1.5	2.0	2.5	3.0	3.5	4.0	4.5	5.0	5.5	6.0	6.5
LUNGS	124	216	301	375	439	493	539	578	610	637	660	678
ACTIVE BONE MARROW	27	49	71	92	113	131	149	165	180	194	207	219
THYROID	4.0	11	21	30	40	49	57	64	70	75	79	82
TRUNK TISSUE	52	82	109	132	152	170	185	199	210	221	230	238
CDI (10 ⁻⁵)	0.65	1.09	1.49	1.85	2.16	2.44	2.68	2.88	3.06	3.22	3.36	3.48
TESTES	+	+	+	+	+	0.1	0.1	0.1	0.1	0.1	0.1	0.2

Organ dose estimation: PCXMC

- A commercial computer program for calculating patients' organ and effective doses in radiography examinations
 - Developed by Tapiovaara et al. (STUK, Finland)
 - Current version, PCXMC 2.0 (released in Nov 2008)
 - Based on the computational phantoms (Cristy and Eckerman, 1987) coupled with Monte Carlo transport technique

Graphical interface for user input of technique factors

Radiography
Mammography
Fluoroscopy
CT

Defform [C:\Program Files\PCXMC\MCRUNS\15-year-old\15y-Abdomen-AP.DF2]

File

Main menu New Form Open Form Save Form Save Form As ... Print As Text

Header text Typical abdomen AP, 15y

Phantom data

Age: ☐ 0 ☐ 1 ☐ 5 ☐ 10 ☒ 15 ☐ Adult

Phantom height 168.10 Phantom mass 55.90 ☒ Arms in phantom

Standard: 168.1 Standard: 56.3

Geometry data for the x-ray beam

FSD 85.40 Beam width 23.32 Beam height 31.04 Xref 0.0000 Yref 0.0000 Zref 19.9900

Projection angle 270.00 Cranio-caudal angle 0.00

LATR=180 AP=270 (pos) Cranial X-ray tube
LATL=0 PA=90 (neg) Caudal X-ray tube

☒ Draw x-ray field

Draw Update Field Stop

MonteCarlo simulation parameters

Max energy (keV) 150 Number of photons 50000

Field size calculator

FID 110 Image width 18 Image height 24 Calculate

Phantom exit- image distance: 5.0

FSD Beam width Beam height Use this data

☒ Skeleton ☒ Brain ☒ Heart ☒ Testes ☒ Spleen ☒ Lungs ☒ Ovaries ☒ Kidneys ☒ Thymus ☒ Stomach ☒ Salivary glands ☒ Oral mucosa

☒ Pancreas ☒ Uterus ☒ Liver ☒ Upper large intestine ☒ Lower large intestine ☒ Small intestine ☒ Thyroid ☒ Urinary bladder ☒ Gall bladder ☒ Oesophagus ☒ Prostate ☒ Pharynx/trachea/sinus

Rotation increment 30 View angle 330

Quick Sharp

Input measurement and organ dose output

Radiography
Mammography
Fluoroscopy
CT

Patient input dose

Input dose value: mGy

Incident air kerma value used in calculations: mGy

[Corresponds to about 9.7mAs]

Input dose quantity and unit:

- ☒ Incident air kerma (mGy)
- ☐ Dose-Area Product (mGycm²)
- ☐ Entrance exposure (mR)
- ☐ Exposure -Area Product (Rcm²)
- ☐ Current -Time Product (mAs)

(Input dose quantities are for measurements without BSF)

OK ! Cancel

PCXMC- Dose Calculation

File Run

Main menu Change X-ray Spectrum Open MC data for dose calculation Print Save As ...

X-ray tube potential: 120 kV Filtration: 3 mm Al + 0,1 mm Cu
Anode angle: 17 deg

File: C:\Program Files\PCXMC\MCRUNS\Adult\adult-chest-PA200000.en2
typical chest PA, Adult Phantom: Adult, Arms included. Simulation: Photons/Energy level: 200000 Maximum energy: 150 keV
Projection angle (LATL=0,PA=90,LATR=180,AP=270): 90.000 Obl. angle: 0.000
Field width: 30.10 cm and height: 37.60 cm FSD: 160.000 cm Ref.point (x,y,z(cm)): (0.000, 0.000, 52.000)
Phantom height: 178.600 cm and mass: 73.200 kg Scaling factors: sx(=sy): 1.000 and sz: 1.000
Incident air kerma:..... 0.150 mGy Tube voltage: 120 kV Filter:..... 3 mm Al + 0.1 mm Cu

Organs	Dose (mGy)	Error (%)	Organs	Dose (mGy)	Error (%)
Active bone marrow	0,052206	0,1	(Scapulae)	0,400557	0,4
Adrenals	0,154045	1,7	(Clavicles)	0,060160	1,4
Brain	0,000822	3,3	(Ribs)	0,288670	0,2
Breasts	0,035366	0,9	(Upper arm bones)	0,035092	1,0
Colon (Large intestine)	0,003454	1,5	(Middle arm bones)	0,039175	1,0
(Upper large intestine)	0,005363	1,6	(Lower arm bones)	0,006518	1,8
(Lower large intestine)	0,000929	4,3	(Pelvis)	0,002112	2,1
Extrathoracic airways	0,007761	4,7	(Upper leg bones)	0,000026	18,1
Gall bladder	0,027439	2,1	(Middle leg bones)	0,000001	47,8
Heart	0,058200	0,6	(Lower leg bones)	0,000000	NA
Kidneys	0,092431	0,6	Skin	0,024164	0,2
Liver	0,065439	0,3	Small intestine	0,004068	1,4
Lungs	0,131340	0,2	Spleen	0,124134	0,8
Lymph nodes	0,035024	0,4	Stomach	0,037965	1,1
Muscle	0,027625	0,1	Testicles	0,000001	69,9
Esophagus	0,081588	1,0	Thymus	0,029615	3,1
Oral mucosa	0,002377	6,2	Thyroid	0,023657	3,7
Ovaries	0,000911	21,3	Urinary bladder	0,000195	17,7
Pancreas	0,073502	1,1	Uterus	0,000731	8,1
Prostate	0,000068	83,9			
Salivary glands	0,004203	3,5	Average dose in total body	0,035716	0,1
Skeleton	0,074937	0,1	Effective dose ICRP60 (mSv)	0,042921	0,3
(Skull)	0,004391	1,5	Effective dose ICRP103 (mSv)	0,045029	0,2
(Upper Spine)	0,044352	1,2			
(Middle Spine)	0,326392	0,3			
(Lower Spine)	0,075329	0,9	Abs. energy fraction (%)	58,042776	

Risk assessment output

Radiography
Mammography
Fluoroscopy

PCXMC - Risk assessment

File Run About

Main menu Open dose data (and clear old doses)... Add further dose data... Calculate risks Clear doses Print report Save report as..

Age: 47.0 Gender: ☒ Male ☐ Female Statistics: ☐ Euro-American ☐ Asian ☒ Finnish

Active bone marrow (mSv) 0.1068
Breasts (women) (mSv) 0
Colon (mSv) 0.008086
Liver (mSv) 0.06977
Lungs (mSv) 0.2945
Ovaries (women) (mSv) 0
Prostate (men) (mSv) 8.4E-5
Stomach (mSv) 0.2233
Thyroid (mSv) 0.08385
Uterus (women) (mSv) 0
Urinary bladder (mSv) 0.000563
Weighted remainder (mSv) 0.1045

Dose files summed:
Sample-chestPA.mG2
Sample-chestLAT.mG2

Input data
Finnish mortality data
47.0 year-old male
Sum of incident air kermas in the selected dose files: 0.58 mGy
Sum of effective doses in the selected dose files : 0.1408 mSv


Stochastic radiation risks
Risk of exposure-induced cancer death (REID): 0.000446 %
(Cancer mortality for other causes; not related to this exposure: 20.1 %)
Expected length of remaining life 28.2 years
Loss of life expectancy (LLE): 0.7 hours
LLE/REID : 17.1 years

Risk of exposure-induced death (REID) for various cancers

Cancer type	REID
leukemia	5.31E-5 %
colon cancer	2.79E-6 %
liver cancer	9.91E-6 %
lung cancer	0.000239 %
prostate cancer	5.7E-9 %
stomach cancer	4.86E-5 %
thyroid cancer	1.73E-6 %
bladder cancer	9.93E-8 %
other cancer	9.05E-5 %

Mammography

- Average (or mean) glandular dose (AGD)
 - Used to describe the dose to the breast
 - Considered to be at greatest risk
 - Replaced traditional quantities (skin dose, midplane breast dose, and etc.)



Factors affecting dose in mammography

- Beam energy
 - Approximately 24-30 kVp
 - Small difference in beam energy affect breast dose
 - Higher beam energy reduce breast dose
- Target material
 - Molybdenum (18 and 20 keV) and rhodium (20 and 23 keV)
 - Rhodium used for thicker breast
- Filter material
 - Molybdenum and rhodium

Factors affecting dose in mammography

- Grids
 - Reduce the scattered radiation to increase image contrast
 - High contrast images are very important because of similar composition of glandular tissue with surrounding ones
- Magnification
 - Move breast closer to the x-ray tube
 - 1.5 to 2.0 times magnified
 - Increase breast dose according to the inverse square law

Factors affecting dose in mammography

- Breast thickness and tissue composition
 - Thick (or large) breasts or those with dense composition need higher energy beam and longer exposure time, and receive higher AGD
- Compression
 - Provides better imaging geometry
 - Lower AGD to the patient
 - More uniform exposure the breast

Organ dose estimation: Conversion factor

- The average glandular dose, $D_g = D_{gN} \times X_{ESE}$

X_{ESE} : the entrance skin exposure (measurable)

D_{gN} : ESE-to-AGD conversion factor (obtained from Monte Carlo simulation)

TABLE 8-6. D_{gN} CONVERSION FACTOR (mRAD PER ROENTGEN) AS A FUNCTION OF HVL AND kVp FOR Mo TARGET/FILTER: 4.5-CM BREAST THICKNESS OF 50% GLANDULAR AND 50% ADIPOSE BREAST TISSUE COMPOSITION*

HVL (mm)	kVp							
	25	26	27	28	29	30	31	32
0.25	122							
0.26	126	128						
0.27	130	132	134					
0.28	134	136	138	139				
0.29	139	141	142	143	144			
0.30	143	145	146	147	148	149		
0.31	147	149	150	151	152	153	154	
0.32	151	153	154	155	156	158	159	160
0.33	155	157	158	159	160	162	163	164
0.34	160	161	162	163	164	166	167	168
0.35	164	166	167	168	169	170	171	172
0.36	168	170	171	172	173	174	175	176
0.37		174	175	176	177	178	178	179
0.38			179	180	181	182	182	183
0.39				184	185	186	186	187
0.40					189	190	191	192

*Adapted from ACR QC Manual, 1999.

Advanced conversion factor*

Table A3. Dose conversion coefficients (this work) derived for typical protocols by time period and compressed breast thickness (CBT).


Period	Technique	Target-filter	HVL ^a or filtration	CBT (cm)	kV		Dose coefficient D _{BN}	
					min	max	min	max
1960-1964	Egan	W-Al	0.9 mm Al inherent	3	22	24	0.305	0.353
				5	26	35	0.258	0.377
				8	26	35	0.169	0.255
	<u>Gershon-Cohen</u>	W-Al	1.5 mm Al 1mm inherent	3	25	30	0.449	0.535
				5	25	30	0.294	0.365
				8	25	30	0.192	0.243
1965-1969	Egan	Mo-Mo	HVL = 0.4	3	26	30	0.305	0.309
				5	26	30	0.190	0.194
				8	26	30	0.123	0.126
			HVL = 0.61	3	26	30	0.449	0.453
				5	26	30	0.290	0.297
				8	26	30	0.189	0.195
			0.78 mm Al	3	26	30	0.347	0.375
				5	26	30	0.217	0.238
				8	26	30	0.141	0.156

* Thierry-Chef et al. (in review)

Period	Technique	Target-filter	HVL ^a or filtration	CBT (cm)	kV		Dose coefficient <u>D_{EN}</u>	
					min	max	min	max
1980-1984	<u>Xeroradiography</u>	W-Al	HVL = 1	3	40	55	0.652	0.643
				5	40	55	0.469	0.470
				8	40	55	0.323	0.329
	Screen Film (Low-Dose)	Mo-Mo	HVL = 0.31	3	28	-	0.252	-
				5	28	-	0.156	-
				8	28	-	0.101	-
1985-1989	<u>Xeroradiography</u>	W-Al	HVL = 1.26	3	44	45	0.714	0.713
				5	44	45	0.524	0.524
				8	44	45	0.365	0.366
	Screen Film (Low-Dose)	Mo-Mo	HVL = 0.37	3	27	29	0.286	0.289
				5	27	29	0.178	0.180
				8	27	29	0.115	0.117
		Mo-Mo	HVL = 0.49	3	27	29	0.368	0.370
				5	27	29	0.232	0.234
				8	27	29	0.151	0.152
1990-1999	<u>Xeroradiography</u>	W-Al	HVL = 1.3	3	46	-	0.726	-
				5	46	-	0.536	-
				8	46	-	0.375	-
	Screen Film (Low-Dose)	Mo-Mo	HVL = 0.35	3	25	28	0.269	0.274
				5	25	28	0.166	0.170
				8	25	28	0.108	0.111
			HVL = 0.37	3	25	28	0.283	0.287
				5	25	28	0.175	0.179
				8	25	28	0.113	0.116
2000+	Screen Film (Low-Dose)	Mo-Mo	0.03 mm Mo	3	24	28	0.241	0.252
				5	24	28	0.149	0.156
				8	24	28	0.097	0.101

□

Fluoroscopy




Factors affecting dose in fluoroscopy

- Beam energy
 - Higher kVp results in more penetrating beam and reduces tube current
- Collimation
 - Use the smallest field to image only the area of interest
 - Reduce the scattered radiation and leads to higher-quality images

Factors affecting dose in fluoroscopy

- Increase source-to-skin distance
 - Reduce the patient dose according to inverse square law
- Decrease patient-to-image intensifier distance
 - Reduce the patient dose since lower x-ray fluence is needed for acceptable image quality
 - Low image quality due to the increased scattered radiation



Factors affecting dose in fluoroscopy

- Image magnification
 - Move the image intensifier farther from the patient or
 - Move x-ray source closer to patient
 - Increase the patient dose
- Grids
 - Reduce the scattered radiation to increase image contrast
 - Patient doses increase by a factor of two or more
- Patient size
 - kVp and tube current must be increased for thicker patients

Factors affecting dose in fluoroscopy

- Beam-on time
 - Directly proportional to the patient dose
 - Several techniques to reduce beam-on time
 - Being aware of the amount of the beam-on time
 - Last-frame-hold feature (display the last image after the beam is off)
 - Aggressive use of low frame rate pulsed fluoroscopy
 - Release the fluoroscopy pedal frequently

Organ dose estimation: conversion factor*

- Heavily rely on computer simulation using Monte Carlo transport technique and computational human phantoms




Table 2. Organ dose conversion coefficients (mGy per Gy cm²) and σ

AP projection	60 kVp 3.5 mm Al UFHADM weight percentile		
	10%	50%	90%
Organs			
Colon	6.04E-03 (0.67%)	5.51E-03 (0.63%)	5.61E-03 (0.63%)
Lung	4.70E-01 (0.05%)	2.83E-01 (0.06%)	1.45E-01 (0.08%)
Stomach	2.33E-01 (0.15%)	1.85E-01 (0.15%)	1.10E-01 (0.19%)
Bladder	8.88E-05 (10.14%)	7.53E-05 (9.21%)	6.52E-05 (10.62%)
Liver	1.74E-01 (0.10%)	1.32E-01 (0.10%)	8.28E-02 (0.13%)
Esophagus	2.56E-01 (0.21%)	1.77E-01 (0.23%)	8.71E-02 (0.32%)
Thyroid	2.34E-02 (1.26%)	2.04E-02 (1.19%)	1.81E-02 (1.25%)
Gonads	4.71E-05 (20.47%)	7.74E-05 (18.29%)	3.31E-05 (24.10%)
Skin	8.07E-02 (0.03%)	7.26E-02 (0.03%)	6.56E-02 (0.03%)
Brain	2.68E-04 (2.77%)	2.64E-04 (2.40%)	2.22E-04 (2.65%)
Kidneys	1.27E-02 (0.62%)	1.00E-02 (0.63%)	6.58E-03 (0.77%)
Salivary glands	2.55E-03 (1.94%)	3.45E-03 (1.46%)	2.40E-03 (1.81%)
Adrenals	4.52E-02 (0.93%)	3.31E-02 (0.96%)	1.95E-02 (1.24%)
Gall bladder	2.90E-02 (0.83%)	2.46E-02 (0.80%)	1.88E-02 (0.91%)


* Johnson et al. PMB (2009)

Dose estimation: Skin dose

- Direct dose measurement
 - Thermoluminescent dosimeter (TLD)
 - X-ray film
- Real-time direct dose measurement
 - MOSFET dosimeter
- Indirect dose measurement
 - Measure dose at the collimator port
 - Dose derived from system parameters (e.g. PEMNET system)
- Real-time parameters
 - Fluoroscopic time
 - Dose-area-product

Radiography
Mammography
Fluoroscopy
CT


Dose estimation: Operator*



* Kim et al. HP (2008)


Radiography
Mammography
Fluoroscopy
CT

Dose estimation: Operator*




* Kim et al. HP (2008)

Dose estimation: Operator*




* Kim et al. HP (2008)


Computed Tomography




Two innovations in CT



Helical scan:
Faster scan time



Multi-detector:
More information




Measurable quantities in CT

- Computed Tomography Dose Index (CTDI)₁₀₀
 - Single axial rotation
 - 100-mm long ion chamber and head/body CTDI phantoms




- Weighted CTDI: $CTDI_w = 1/3 CTDI_{100,center} + 2/3 CTDI_{100,peripheral}$
- Volume-weighted CTDI: $CTDI_{vol} = CTDI_w / pitch$
- Dose Length Product (DLP) = $CTDI_{vol} \times \text{scan length (cm)}$

- Not designed for or representing patient organ dose!




Pitch in helical scan



Volume-weighted CTDI: $\text{CTDI}_{\text{vol}} = \text{CTDI}_w / \text{pitch}$

Factors affecting dose in CT

- Beam energy
- Photon fluence (current-time-product)
- Helical pitch
- Patient size



CTDI body (left) and head (right) phantoms

Factors affecting dose in CT: Energy

Table 1
Changes in CTDI_w in Head and Body
Phantoms as a Function of Kilovolt Peak

Beam Energy (kVp)	CTDI _w in Head Phantom (mGy)	CTDI _w in Body Phantom (mGy)
80	14	5.8
100	26	11
120	40	18
140	55	25

Note.—All other factors were held constant at 300 mA, 1 sec, and 10 mm. Results are from a single-detector CT scanner.

$$14 \times \left(\frac{140}{80} \right)^{2.5} = 56.7$$

Factors affecting dose in CT: Fluence (mAs)

Table 2
Changes in CTDI_w in Head and Body
Phantoms as a Function of Milliampere-
Seconds Setting

Tube Current– Time Product (mAs)	CTDI _w in Head Phantom (mGy)	CTDI _w in Body Phantom (mGy)
100	13	5.7
200	26	12
300	40	18
400	53	23

Note.—All other factors were held constant at 120 kVp and 10 mm. Results are from a single-detector CT scanner.


Factors affecting dose in CT: Pitch

Table 3
Changes in CTDI_{vol} in Head and Body
Phantoms as a Function of Pitch


Pitch	CTDI _{vol} in Head Phantom (mGy)	CTDI _{vol} in Body Phantom (mGy)
0.5	80	36
0.75	53	24
1.0	40	18
1.5	27	12
2.0	20	9

Note.—All other factors were held constant at 120 kVp, 300 mA, 1 sec, and 10 mm. Results are from a single-detector CT scanner.

Factors affecting dose in CT: Patient size*




16-cm diameter head phantom



32-cm diameter head phantom

CTDI_{vol} vs. actual organ dose

Organ dose per CTDI_{vol} (abdomen-pelvis scan for adult male)*




* Lee et al. Medical Physics (2011)

Organ dose estimation: Software tools

ImPACT

- NRPB database (UK)
- Hermaphrodite adult
- **No children**



ORNL adult
hermaphrodite
phantom

CT-Expo


- GSF database (Germany)
- Male and female adult
- **Two children**



BABY



CHILD



ADAM

EVA

ImPACT

ImPACT CT Patient Dosimetry Calculator

Version 1.0.2 12/11/2009

Scanner Model:

Manufacturer: Siemens

Scanner: Siemens Sensation 16

kV: 120

Scan Region: Head

Data Set: MCSET21

Current Data: MCSET21

Scan range

Start Position: 0 cm

End Position: 43 cm

Acquisition Parameters:

Tube current: 100 mA

Rotation time: 1 s

Spiral pitch: 1

mAs / Rotation: 100 mAs

Effective mAs: 100 mAs

Collimation: mm

Rel. CTDI: Look up 1.00 (assumed)

CTDI (air): Look up 21.8 mGy/100mAs

CTDI (soft tissue): 23.3 mGy/100mAs

n CTDI_w: Look up 16.6 mGy/100mAs

Organ weighting scheme: ICRP 60

Organ	w_T	H_T (mGy)	$w_T \cdot H_T$
Gonads	0.2	6.5	1.3
Bone Marrow	0.12	4.8	0.57
Colon	0.12	11	1.3
Lung	0.12	1.2	0.15
Stomach	0.12	11	1.3
Bladder	0.05	12	0.61
Breast	0.05	0.35	0.018
Liver	0.05	10	0.5
Oesophagus (Thymus)	0.05	0.26	0.013
Thyroid	0.05	0.031	0.0016
Skin	0.01	3.3	0.033
Bone Surface	0.01	5.8	0.058
Not Applicable	0	0	0
Not Applicable	0	0	0
Remainder	0.025	4.7	0.12
Kidneys	0.025	13	0.31
Total Effective Dose (mSv)			6.3


Remainder Organs	H_T (mGy)
Adrenals	9.4
Small Intestine	11
Kidney	13
Pancreas	9.4
Spleen	10
Thymus	0.26
Uterus	12
Muscle	4.7
Brain	0.0016
Not Applicable	N/A
Not Applicable	N/A
Not Applicable	N/A
Not Applicable	N/A
Other organs of interest	H_T (mGy)
Eye lenses	0.0027
Testes	2.3
Ovaries	11
Uterus	12
Prostate	12

CTDI_w: 16.6 mGy

CTDI_{vol}: 16.6 mGy

DLP: 713 mGy.cm


Scan Description / Comments



CT-EXPO

Calculate

2. Scan Range



Get Values

Scan Range Data (Slice Positions)

Scan Range z from z-	to z+	L [cm]
40	70	30

1. Age Group Gender

Adult

☒ male
 ☐ female

3. Scanner Model

Manufacturer

General Electric

Scanner

9800

Scanner Data for Scan Region "Body"

$nCTDI_w$ [mGy/mAs]	U_{ref} [kV]	$P_{B,H}$	k_{CT}	k_{OB}	ΔL [cm]
0.062	120	0.24	0.50	1.00	0.0

4. Select mode

☐ Body mode for head/neck region
 ☒ Spiral mode

5. Scan Parameters

Please Enter Actual Settings:

U [kV]	I [mA]	t [s]	Q_{el} [mAs]	Q [mAs]	$N * h_{col}$ [mm]	TF [mm]	h_{rec} [mm]	p	Ser.
120	120	2	240	0	10.0	10.0	1.0	1.0	1

6. Results

Dose Values per Scan or per Series*

$CTDI_w$ [mGy]	$CTDI_{vol}$ [mGy]	DLP_w^* [mGy*cm]	E^* [mSv]	D_{uterus}^* [mSv]
14.9	14.9	448	6.4	n.a.

Child/Baby: all CTDI and DLP values refer to 16cm head phantom !

Dose Values per Examination

DLP_w [mGy*cm]	E [mSv]	D_{uterus} [mSv]
448	6.4	n.a.


Effective dose E refers to ICRP 60

Please note:


All organ doses H_T are based on conversion coefficients for stand-ard patients (ADAM, EVA, CHILD, BABY) and serve for information purposes only (in particular organs outside the scan range) !

Tissue or Organ	H_T per Series [mSv]	Remainder Organs	H_T per Series [mSv]
Thyroid	6.4	Brain	0.3
Breasts	0.0	Thymus	22.7
Oesophagus	22.7	Spleen	5.8
Lungs	23.4	Pancreas	5.2
Liver	7.4	Adrenals	9.3
Stomach	4.8	Kidneys	1.5
Colon	0.0	Small intest.	0.2
Testicles	0.0	Upp. large int.	0.3
Ovaries	0.0	Uterus	0.0
Bladder	0.0		
Bone marrow	5.6	Misc.	H_T per Series [mSv]
Bone surface	15.0		
Skin	5.4	Eye lenses	0.4

CT-EXPO



Organ dose estimation: Perimeter-base*



Perimeter-based organ doses*

TABLE II. Results of exponential regression analysis describing $\overline{nD}_{P,O}$ as a function of perimeter (cm) for fully irradiated organs.

Organs	Exponential regression coefficients		Correlation coefficient
	A_O	B_O	R^2
Liver	3.824	−0.0120	0.98
Stomach	3.780	−0.0113	0.97
Adrenals	4.029	−0.0128	0.95
Kidney	3.969	−0.0124	0.99
Pancreas	3.715	−0.0122	0.97
Spleen	3.514	−0.0111	0.95
Gall bladder	3.994	−0.0115	0.95

Summary

- Epidemiology needs individualized organ dose.
- Three approaches
 - Measurement: expensive, labor-intensive, and not individualized
 - Calculation: cost-effective, fewer man-hour, and individualized
 - Conversion factor: derived from calculation
- Four different imaging modalities
 - Radiography
 - Mammography
 - Fluoroscopy
 - Computed Tomography

References

- Parry RA, Glaze SA, and Archer BR 1999 The AAPM/RSNA Physics Tutorial for Residents. Radiographics **19(5) 1289-1302**
- McNitt-Gray MF 2002 AAPM/RSNA physics tutorial for residents: topics in CT. Radiographics **22(6) 1541**

Thank you for your attention!
Any questions or comments appreciated